


In Vitro (modus 2)

for glass og forhåndsprosessert lydfil

2008

Ragnhild Berstad



NB
noter

this music is copyright protected

IN VITRO (modus 2)

GLASS
ampl
+ rev

LYDFIL

-15" -10" -5" 0" 5" 10" 15" 20" 25" 30"

30" 35" 40" 45" 50" 55" 1:00 1:05 1:10 1:15

1:15 1:20 1:25 1:30 1:35 1:40 1:45 1:50 1:55 2:00

sfz

GLASS
amp
+ reb

LYDFIL

Ragnhild Berstad 5/2 2008

Tuning of crystal glasses

glasses in need to be tuned, if borrowed from composer

tune glasses with water

Glasses are to be amplified

Add a slight reverberation, approx 2.5 sec rev. time